

## LAMPS AND LIGHTS

We service nearly every make and model in the dental care industry!

### LED Blast Fast Curing Light

The LED Blast fast curing light utilizes a chip-on-board LED module to create exceptional performance.

LED technology eliminates costly halogen bulbs and reduces heat buildup. The new LED Blast cures in less than 10 seconds with 1000 to 1300 mw/cm2 of continuous power. Comes complete with 11mm autoclavable light guide, anti-glare shield, and charger base with power supply. Long life batteries allow use all day and provides 750 five-second cures with an overnight charge. Timing cycles of 5, 10, 20 and 30 seconds. Ergonomic pistol shape can be used with either hand.



- LB001** ..... LED Blast
- L009-7000** ..... Hilux Light Meter

### LIGHT GUIDES

- M014-01** ..... 3mm, 60 deg
- M011-24** ..... 8mm, 60 deg
- M015-24** ..... 11mm, 60 deg

### Blast™-Lite and Mini-Blast™



Wand shaped lights such as the Blast™-Lite and Mini-Blast™ are smaller and lighter than most hand-held curing lights, ideal for quick procedures. Both high power lights are equipped with timer settings of 5, 10 and 20 seconds. Included are long lasting Lithium-Ion batteries for extended use. Removable glass tip does not require a sleeve and can be sterilized, metal case allows for easy cleaning.

- LB005** ..... Blast™-Lite
- LB006** ..... Mini-Blast™

phone **800.437.5606** | fax **800.235.9001** | email **sales@nexadental.com**  
order online anytime at **www.nexadental.com**



Orders ship UPS ground. US dollars, checks, Visa, MasterCard and American Express accepted. All orders FOB La Vista, NE. Prices subject to change. Custom products may not be returned for refund or exchange.

## Replacement Bulbs for Operatory, Curing, Fiber Optic, Handpiece, Video

NEXADENTAL has the special lamps you need to keep your equipment operating in top form. Don't see what you need here? Contact your sales representative to find the lamp you need!



<b>BW.BHB</b> .....	120V,250W	<b>BW.13165</b> .....	14V,35W
<b>BW.DNE</b> .....	120V,150W	<b>BW.13165-IMP</b> .....	14V,35W
<b>BW.DNF</b> .....	21V,150W	<b>BW.3367-5267</b> .....	10V,52W
<b>BW.DZE</b> .....	24V,150W	<b>BW.13865</b> .....	12V,75W
<b>BW.EFM</b> .....	8V,50W	<b>BW.21237</b> .....	12V,80W
<b>BW.EFN</b> .....	12V,75W	<b>BW.3367-5267</b> .....	10V,52W
<b>BW.EFP</b> .....	12V,100W	<b>BW.42001021GEN1</b> ....	DABI ATLANTE
<b>BW.EFP-OSR</b> .....	12V,100W	<b>BW.42001021GEN3</b> ....	DABI ATLANTE
<b>BW.EFR</b> .....	15V,150W	<b>BW.553-3881</b> .....	2 Pack Kavo Lux Coupler Bulb
<b>BW.EHJ</b> .....	24V,250W	<b>BW.561398</b> .....	13V,100W
<b>BW.EHT</b> .....	120V,250W	<b>BW.60986</b> .....	2 Pack Star Quick Connect Swivel
<b>BW.EJA</b> .....	21V,150W	<b>BW.64115</b> .....	12V,20W
<b>BW.EJL</b> .....	24V,200W	<b>BW.64440</b> .....	12V,50W
<b>BW.EJM</b> .....	21V,150W	<b>BW.64602</b> .....	12V,50W
<b>BW.EJV</b> .....	21V,150W	<b>BW.64613</b> .....	12V,75W
<b>BW.EKE</b> .....	21V,150W	<b>BW.64624</b> .....	12V,100W
<b>BW.ELZ</b> .....	21V,150W	<b>BW.650064-GEN2</b> .....	2 Pack Midwest Power Optic
<b>BW.ENL</b> .....	12V,50W	<b>BW.650197</b> .....	3 Pack Midwest Power Optic
<b>BW.ENZ</b> .....	30V,50W	<b>BW.70201102004</b> .....	13V,75W
<b>BW.EPV</b> .....	14.5V,90W	<b>BW.95775</b> .....	41V,200W
<b>BW.EPX</b> .....	14.5V,90W	<b>BW.JA12V-55WD</b> .....	12V,55W
<b>BW.ESD</b> .....	120V,150W	<b>BW.JA12V-100WD</b> .....	12V,100W
<b>BW.ESJ</b> .....	82V,85W	<b>BW.JA24V-70WD</b> .....	24V,70W
<b>BW.ESJ-GE</b> .....	82V,85W	<b>BW.JC12V-50W</b> .....	12V,50W
<b>BW.ESP</b> .....	120V,150W	<b>BW.JC12V-50WG1</b> .....	12V,50W
<b>BW.ESX</b> .....	12V,20W	<b>BW.JCD24V-75W</b> .....	24V,75W
<b>BW.ESX-13.8V</b> .....	13.8V,25W	<b>BW.JCD120V-150WB</b> ..	120V,150W
<b>BW.ESY</b> .....	100V,150W	<b>BW.JCR120V150WB</b> ...	120V,150W
<b>BW.ETJ</b> .....	120V,250W	<b>BW.JCRM12V-50W</b> .....	12V,50W
<b>BW.EVA</b> .....	12V,100W	<b>BW.JCRM12V-75W</b> .....	12V,75W
<b>BW.JCR120V150WB</b> ..	120V,150W	<b>BW.JCRM12V-100W</b> ...	12V,100W
<b>BW.FCS</b> .....	24V,150W	<b>BW.M34</b> .....	6V,20W
<b>BW.FDS.DZE</b> .....	24V,150W	<b>BW.Q75CL.DC-28V</b> .....	28V,75W
<b>BW.FDT</b> .....	12V,100W	<b>BW.Q150T4CL25IM</b> ....	25V,150W (Import Brand)
<b>BW.FTD</b> .....	12V,20W	<b>BW.Q150T4CL25GE</b> ....	25V,150W (GE Brand)
<b>BW.25A-120V</b> .....	120V,25W	<b>BW.THE.71</b> .....	15V,150W
<b>BW.100A21SP</b> .....	120V,100W	<b>BW.WA01100</b> .....	2.5V VACUUM LAMP
<b>BW.125T10.P</b> .....	120V,125W	<b>BW.Y900529-GEN1</b> .....	NSK HANDPIECE BULB
<b>BW.300SB.IF</b> .....	120V,300W	<b>BW.Y900529-GEN3</b> .....	NSK HANDPIECE BULB
<b>BW.041-317-00</b> .....	ADEC Handpiece Bulb		

## Replacement Bulbs for X-Ray Viewing; Magnifying Lamps; Work Lights

NEXADENTAL has the special lamps you need to keep your equipment operating in top form. Don't see what you need here? Contact your sales representative to find the lamp you need!



- BW.CF9DS.E.841** ..... 9W Compact Fluorescent 4100K.
- BW.DS9W.78**..... 9W Blue Compact Fluorescent
- BW.F6T5.CW** ..... 6W 9" Cool White
- BW.F6T5.D**..... 6W 9" Daylight
- BW.F8T5.CW** ..... 8W 12" Cool White
- BW.F8T5.D**..... 8W 12" Daylight
- BW.F13T5.CW** ..... 13W 21" Cool White
- BW.F14T8.CW** ..... 14W 15" Cool White 1" Diameter
- BW.F14T8.D**..... 14W 15" Daylight 1" Diameter
- BW.F14T12.CW** ..... 14W 15" Cool White 1 1/2" Diameter
- BW.F14T12.D** ..... 14W 15" Daylight 1 1/2"Diameter
- BW.F15T8.CW** ..... 15W 18" Cool White
- BW.F15T8.D** ..... 15W 18" Daylight
- BW.F30T8.C50** ..... 30W Color Corrected
- BW.FC8T9.CW** ..... 22W 8" Diameter Ring Fluorescent

## Replacement Bulbs for Indicator, Panel, Communicator & Miniature Lamps

NEXADENTAL has the special lamps you need to keep your equipment operating in top form. Don't see what you need here? Contact your sales representative to find the lamp you need!



- BW.47** ..... 6.3V,.15A Miniature Bayonet Base
- BW.85**..... 28V,.04A Wedge Base
- BW.1206PSB5** ..... 120V,.025A Slide Base
- BW.336**..... 14V,.08A Midget Groove
- BW.387** ..... 28V,.04A Midget Groove
- BW.787** ..... 6V,10W Bi Pin Base
- BW.1073** ..... 12.8V,1.80A SC Bayonet Base
- BW.462420** ..... 24V,20MA Slide Lamp
- BW.552450** ..... 24V,50MA Slide Lamp
- BW.7382** ..... 14V,.08A T1 3/4 Bi Pin Base
- BW.HPR40** ..... 6V,0.67A T3 SC Mini Flange
- BW.HPR53** ..... 4V,.85A T3 SC Mini Flange